

REICHERT® TONO-VERA® TONOMETER IS NOW AVAILABLE!



The Tono-Vera® Tonometer features the advanced, patented **ActiView™ Positioning System** that quickly guides you to the apex of the cornea and automatically takes measurements, giving you confidence in your intraocular pressure (IOP) readings.



Quick and Easy to Use

Objective measurements can be achieved in as little as 3 readings.



Secure Communication

Transfer data to your EMR software via Bluetooth, quickly and securely.



Utmost Comfort for Patients

No anesthesia drops are required.



Easy Storage

Comes with a docking base that stores both Tono-Vera and Ocu-Dot probes.

[LEARN MORE](#)