Measurements made easy.





eye-ruler™ 2 is a well-rounded digital measuring application.

With its intuitive navigation, eye-ruler™ 2 is designed to take all necessary measurements in less than one minute.

It provides practitioners with a complete solution to expertly deliver personalized lenses and ensure a seamless optical experience for the consumer.

1. Embrace a digital in-store experience

Seamless navigation

Quickly master fitting parameters and enjoy an easy measurement process thanks to its intuitive interface.


By following the on-screen guidance, the digital display provides reliable results for easy reading and analysis.





For your high-tech practice

The Frame Comparison Module allows the patient to compare and choose their favorite frame by taking pictures of their try-ons to facilitate comparison.



2. Benefit from all necessary measurements

For all far vision measurements


Thanks to 3D reconstruction technology, all standard and personalized measurements* in less than one minute.

Efficient, well-rounded solution

Based on efficient and robust hardware and software elements, eye-ruler™ 2 is a compact and handheld solution.

With dedicated accessories

Measuring clip and its correct positioning on the patient's frame





*Complete list of all parameters available at the end of this brochure.

3. Raise your game with Near Vision Behavior¹

A unique² lens personalization experience

Exclusive³ to Varilux® XR Series™ lenses, the Near Vision Behavior measurement assesses the patient's visual needs and optimize the shape and position of vision volume.





A 4-in-1 in-store measurement

The natural posture and visual behavior in near vision of the patient is captured by recreating a dynamic reading task of daily life.

Varilux® XR Series™ is a tailor-made design from the combination of four parameters for a smooth and complete adaptation to the new lenses.

¹Only available with Varilux* XR Track, XR 4D and XR Pro depending on local market availability.

²Only available to eye care professionals dispensing Essilor* lens innovations compatible with an Essilor Instruments digital device to do the measurements.

³Exclusivity refers to all those products, measurements and technologies covered by patents owned by Essilor*.

Technical specifications

Accessories

Essilor clip



Personalized measurements

Fitting parameters

Pupillary distances & Monocular PDs in Far Vision

Fitting heights in Far Vision

ELD (Eye Lens Distance = vertex distance)

Pantoscopic tilt

Wrap angle

Additional parameters

Reading distance
Effective lens diameter
Segment heights
Near Vision PDs
Head cape
Near Vision Behaviour

 ϵ

eye-ruler2 Deviceless-Brochure-EN-US-Feb2025

ESSILOR INSTRUMENTS USA 8600 W. Catalpa Avenue, Ste. 703 Chicago, IL 60656 USA Made in France.

*As improvements are made, information and data in this product presentation are not contractually-binding and may be modified without prior notice. Please check the latest version with your local Essilor Instruments representative.

www.EssilorInstrumentsUSA.com 855.393.4647 info@essilorinstrumentsusa.com © Essilor International - February 2025

